



WellStart II BST DNA polymerase

Product Number: BS02

Shipping and Storage

-20°C.

Components

Component	BS02
	1600U
WellStart II BST DNA polymerase(8U/μL)	200μL
10×WellStart II BST Reaction Buffer	1.5mL
100mM MgSO ₄ Solution	1.5mL

Description

WellStart II BST DNA polymerase is a recombinant enzyme expressed and purified by Escherichia coli. Its gene is derived from Bacillus stearothermophilus and has undergone partial point mutations based on the original sequence. This protein has stronger 5'→3' DNA polymerase activity, chain displacement activity, reverse transcriptase activity, and no 5'→3' exonuclease activity. Applied to DNA/RNA isothermal amplification (LAMP), multiple displacement amplification (MDA), whole genome amplification (WGA), etc.

Unit definition

The amount of enzyme required to add 10nmol of deoxynucleotides to acidic insoluble substances within 30 minutes at 65°C is defined as 1 active unit (U).

Heat Inactivation

Incubate at 80°C for 5 minutes before inactivation.

Quality Control

After multiple column purification, the purity was detected by SDS-PAGE to be greater than 98%; No exogenous nuclease activity was detected.

Protocol

Isothermal Amplification (LAMP) Operation Guidelines:

Mix the following components in proportion and incubate at 65°C for 30-60 minutes. Incubate at 80°C for 5 minutes to inactivate.

Components	25μL reaction system	Final Concentration
10×WellStart II BST Reaction Buffer	2.5 μL	1× (contain 2mM MgSO ₄)
100mM MgSO ₄ Solution	1.5 μL	6 mM(8mM in total)
dNTP Mix (10mM)	3.5 μL	1.4 mM each
FIP/BIP Primers (25×)	1 μL	1.6 μM
F3/B3 Primers (25×)	1 μL	0.2 μM
LoopF/B Primers (25×)	1 μL	0.4 μM
WellStart II BST DNA polymerase(8U/μL)	0.5-1 μL	160-320 U/mL
DNA Sample	variable	> 10 copies or more
Sterile water	Add to 25 μL	



Tinzyme Co., Limited

Email: sales@tinzyme.com

Website: www.tinzyme.com

Tel: +86-755-86134126

WhatsApp/Facebook/Twitter: +86-189-22896756

Total Volume

25 μ L

Note:1)LAMP primers consist of 4 or 6(including Loop) primers, 25 \times Primers include:40 μ M FIP, 40 μ M BIP, 5 μ M F3, 5 μ M B3,

10 μ M LoopF, 10 μ M LoopB;

2)To optimize the reaction, the Mg²⁺ concentration (4-10mM), enzyme quantity, or reaction temperature (60-72°C) can be adjusted. The optimal reaction temperature for this enzyme is 65-68°C;

3)Do not shake vigorously. Mixing vigorously can inactivate the enzyme;

4)Ensure that there are no bubbles in the reaction system after adding the system.