



Bst 3.0 Enzyme Mix

Product Number: BSM03

Shipping and Storage

-20°C.

Components

Component	BSM03
	200µL
Bst 3.0 Enzyme Mix	200µL
10×Bst 3.0 Reaction Buffer	1.5mL
100mM MgSO ₄ Solution	1.5mL

Description

Bst 3.0 Enzyme Mix is a mixed enzyme containing Bst DNA polymerase and high-temperature resistant reverse transcriptase. Bst DNA polymerase is a recombinant enzyme expressed and purified through E. coli, which undergoes partial point mutations on the basis of the original sequence. It has stronger 5'→3' DNA polymerase activity, strand displacement activity, reverse transcriptase activity, and no 5'→3' exonuclease activity. High temperature resistant reverse transcriptase is a new enzyme modified by genetic engineering, with fast cDNA synthesis speed and significantly improved thermal stability. It can withstand reaction temperatures up to 60°C and is suitable for reverse transcription reactions of RNA templates with complex secondary structures. Bst 3.0 Enzyme Mix can be applied to isothermal amplification reactions (LAMP/RT-LAMP) using RNA or DNA as templates. Bst 3.0 Enzyme Mix can be applied to isothermal amplification reactions (LAMP/RT-LAMP) using RNA or DNA as templates.

Unit definition

This product is suitable for various isothermal amplification reactions such as RT LAMP, LAMP, RCA, CPA, etc.

Heat Inactivation

Incubate at 80°C for 5 minutes before inactivation.

Protocol

Guidelines for Isothermal Amplification (LAMP/RT LAMP) Operation:

Mix the following components in proportion and incubate at 60°C for 30-60 minutes. Incubate at 80°C for 5 minutes to inactivate.

Components	25µL reaction system	Final Concentration
10×Bst 3.0 Reaction Buffer	2.5 µL	1×(contain 2mM MgSO ₄)
100mM MgSO ₄ Solution	1.5 µL	6 mM(8mM in total)
dNTP Mix (10mM)	3.5 µL	1.4 mM each
Primer Mix (25×)	1 µL	
Bst 3.0 Enzyme Mix	0.5-1 µL	
DNA /RNA Sample	variable	
DNA Sample	Add to 25 µL	
Sterile water	25 µL	

Note:1) LAMP primers consist of 4 or 6 (including Loop) primers, 25×Primers include: 40µM FIP, 40µM BIP, 5µM F3, 5µM B3, 10µM LoopF, 10µM LoopB;

2) To optimize the reaction, the Mg²⁺ concentration (4-10mM), enzyme quantity (0.25-1.5µL), or primer concentration can

For Research Use Only



Tinzyme Co., Limited

Email: sales@tinzyme.com

Website: www.tinzyme.com

Tel: +86-755-86134126

WhatsApp/Facebook/Twitter: +86-189-22896756

be adjusted;

3) Do not shake vigorously. Mixing vigorously can inactivate the enzyme;

4) Ensure that there are no bubbles in the reaction system after adding the system.