



WellStart Bst 3.0 Enzyme Mix

Product Number: BSM03W

Shipping and Storage

-20°C.

Components

Component	BSM03W
WellStart Bst 3.0 Enzyme Mix	200µL
10×Bst 3.0 Reaction Buffer	200µL
100mM MgSO ₄ Solution	1.5mL
	1.5mL

Description

WellStart Bst 3.0 Enzyme Mix is a mixed enzyme that contains WellStart Bst DNA polymerase and high-temperature resistant reverse transcriptase. WellStart Bst 3.0 Enzyme Mix adds modifications on top of Bst 3.0, eliminating non-specific amplification generated during reaction establishment at room temperature and eliminating the need for a separate activation step, which can improve the specificity of the reaction. High temperature resistant reverse transcriptase is a new enzyme modified by genetic engineering, with fast cDNA synthesis speed and significantly improved thermal stability. It can withstand reaction temperatures up to 60°C and is suitable for reverse transcription reactions of RNA templates with complex secondary structures. WellStart Bst 3.0 Enzyme Mix can be applied to isothermal amplification reactions (LAMP/RT LAMP) using RNA or DNA as templates.

Application

This product is suitable for various isothermal amplification reactions such as RT LAMP, LAMP, RCA, CPA, etc.

Thermal inactivation

Incubate at 80°C for 5 minutes before inactivation.

Protocol

Guidelines for Isothermal Amplification (LAMP/RT LAMP) Operation:

Mix the following components in proportion and incubate at 60°C for 30-60 minutes. Incubate at 80°C for 5 minutes to inactivate.

Component	25µL reaction system	Final Concentration
10×WellStart Bst 3.0 Enzyme Mix Buffer	2.5µL	1×(contain 2mM MgSO ₄)
100 mM MgSO ₄ Solution	1.5µL	6mM(8mM in total)
dNTP Mix (10mM)	3.5µL	1.4mM each
Primer Mix (25×)	1µL	
WellStart Bst 3.0 Enzyme Mix	0.5-1µL	
DNA /RNA Sample	variable	
Sterile water	Supplemented to 25µL	
Total Volume	25µL	

Note:1)Primers consisting of 4 or 6 (including Loop) primers, 25 × Primers include:40µM FIP, 40µM BIP, 5µM F3, 5µM B3, 10µM LoopF, 10µM LoopB;

2)To optimize the reaction, the Mg²⁺ concentration (4-10mM), enzyme quantity (0.25-1.5µL), or primer concentration can be adjusted;



Tinzyme Co., Limited

Email: sales@tinzyme.com

Website: www.tinzyme.com

Tel: +86-755-86134126

WhatsApp/Facebook/Twitter: +86-189-22896756

- 3) Do not shake vigorously. Mixing vigorously can cause enzyme inactivation;
- 4) Ensure that there are no bubbles in the reaction system after adding the system.