

## Enterovirus 71/Coxsackievirus A16 (EV71CA16) Nucleic Acid

### Detection Kit (Dual Fluorescent PCR Method)

Product Number:DTK515

---

#### Shipping and Storage

1.  $-20^{\circ}\text{C} \pm 5^{\circ}\text{C}$ , stored in the dark, transported, and subjected to repeated freeze-thaw cycles no more than 5 times, with a validity period of 12 months.
2. The above specimens can be stored at  $-20^{\circ}\text{C}$  in the short term and  $-70^{\circ}\text{C}$  in the long term, but cannot exceed 6 months. The specimens should be transported in  $2-8^{\circ}\text{C}$  ice packs, and repeated freezing and thawing are strictly prohibited.

#### Component

Component	50T
EV71/CA16 reaction solution	500 $\mu\text{L} \times 2$
Enzyme solution	50 $\mu\text{L}$
EV71/CA16 positive quality control product	250 $\mu\text{L}$
Negative quality control product	250 $\mu\text{L}$

**Note: Different batches of reagents cannot be mixed.**

#### Description

Hand, foot, and mouth disease is a common childhood illness caused by enteroviruses, mostly resulting in fever and water rash on the mouth, palms, and soles of the feet. A small number of patients may experience symptoms such as meningitis, encephalitis, neurogenic pulmonary edema, myocarditis, and even death. The enteroviruses that cause hand, foot, and mouth disease include Coxsackievirus, novel enteroviruses, and enteroviruses. Currently, laboratory diagnostic methods mainly include virus isolation and nucleic acid testing. Hand foot mouth disease is an acute infectious disease caused by enteroviruses (Coxsackie A group 16, enterovirus 71 is the most common). It mostly occurs in preschool children, especially in the age group under 3 years old, with the highest incidence rate.

This kit is suitable for detecting enterovirus EV71 and Coxsackievirus CA16 RNA in stool or diarrhea samples, and is suitable for auxiliary diagnosis of enterovirus EV71 and Coxsackievirus CA16 infections. The test results are for reference only.

#### Application

This kit uses TaqMan probe method for real-time fluorescence PCR technology, designs a pair of enterovirus EV71 and Coxsackievirus CA16 specific primers, and combines them with a specific probe to amplify and detect the RNA of enterovirus EV71 and Coxsackievirus CA16 in vitro using fluorescence PCR technology, which is used for pathogen diagnosis of suspected infected individuals in clinical practice.

#### Applicable instruments

ABI7500, Agilent MX3000P/3005P, LightCycler, Bio Rad, Eppendorf and other series of fluorescence quantitative PCR detectors.

#### Specimen collection

Symptoms of upper respiratory tract infection (fever, spray, runny nose, cough, etc.): Dip a sampling swab into nasal or pharyngeal secretions and place them in a sterile sampling tube (containing 0.4ml sterile physiological saline)

Fecal: Dip a sampling swab into fecal or diarrheal material and place it in a bacterial sampling tube (containing 0.4ml of sterile physiological saline)

## Protocol

### 1. Sample processing (sample processing area)

#### 1.1. Sample Preparation

Shake the throat swab thoroughly and mix well to obtain the supernatant as the nucleic acid detection material. Take soybean particle sized feces or 100uL water sample abdomen and place it in a 1.5mL centrifuge tube. Extract according to the instructions of the nucleic acid extraction kit.

#### 1.2. Nucleic acid extraction

We recommend using our nucleic acid extraction or purification reagents (magnetic bead method or centrifugal column method) for nucleic acid extraction. Please follow the instructions in the reagent manual.

### 2. Reagent preparation (reagent preparation area)

Based on the total number of samples to be tested, the required number of PCR reaction tubes is N (N=number of samples+1 negative control tube+1 positive control tube); For every 10 portions of reaction tubes, an additional 1 portion is prepared. The preparation of each test reaction system is shown in the following table:

reagent	EV71/CA16 Reaction solution	Enzyme solution
Dosage (sample size N)	19μL	1μL

Transfer the mixed test reaction solution into a PCR reaction tube at a concentration of 20uL per tube.

### 3. Sample addition (sample processing area)

Take 5μL of the nucleic acid, positive control sample, and negative control sample extracted in step 1, and add them to the corresponding reaction tubes. Cover the tubes, mix well, and briefly centrifuge.

### 4. PCR amplification (nucleic acid amplification zone)

4.1. Place the reaction tube to be tested in the reaction tank of the fluorescence quantitative PCR instrument;

4.2. Set the channel and sample information, and set the reaction system to 25μL;

Fluorescence channel selection: Detection channel (Reporter Dye) FAM, Quencher Dye NONE, please do not select ROX reference fluorescence for ABI series instruments, select None.

4.3. Recommended loop parameter settings:

step	Cycles	Temperature	Time	Collect fluorescence signals
1	1 cycle	50°C	10min	No
2	1 cycle	95°C	2min	No
3	45 cycles	95°C	15sec	No
		60°C	30sec	Yes

### 5. Result analysis and judgment

#### 5.1. Result Analysis Condition Setting

(Please refer to the user manuals of each instrument for setting up, taking the ABI7500 instrument as an example)

After the reaction is complete, the results will be automatically saved. Based on the analyzed image, adjust the Start value, End value, and Threshold value of the baseline (users can adjust them according to their actual situation, with Start value set between 3-15 and End value set between 5-20, so that the threshold line is in the exponential period of the amplification curve, and the amplification curve of negative quality control products is flat or below the threshold line). Click Analyze to automatically obtain the analysis results.

#### 5.2. Result judgment

FAM channel shows the detection results of enterovirus EV71; The VIC channel shows the results of Coxsackievirus CA16 detection.

Positive: The Ct value of the detection channel is  $\leq 40$ , and the curve shows a significant exponential growth curve;

Negative: The sample test result shows no Ct value and no specific amplification curve.

Suspicious: If the sample test result is  $40 < \text{Ct value} \leq 45$ , it is recommended to repeat the test. If the detection channel is still  $40 < \text{Ct value} \leq 45$  and the curve has a clear exponential growth curve, it is judged as positive. Otherwise, it is judged as negative.

**Quality control standards**

Negative quality control product: no specific amplification curve or Ct value display;

Positive quality control product: The amplification curve shows a significant exponential growth period, and the Ct value is  $\leq 32$ ;

The above conditions should be met simultaneously, otherwise the experiment will be considered invalid.

**Limitations of detection methods**

1. The results of sample testing are related to the quality of sample collection, processing, transportation, and preservation;
2. Failure to control cross contamination during sample extraction can result in false positive results;
3. Leakage of positive controls and amplification products can lead to false positive results;
4. During the epidemic, genetic mutations and recombination of pathogens can lead to false negative results;
5. Different extraction methods have differences in extraction efficiency, which can lead to false negative results;
6. Improper transportation, storage, or inaccurate preparation of reagents can lead to a decrease in reagent detection efficiency, resulting in false negatives or inaccurate quantitative testing results;
7. The test results are for reference only. If a diagnosis is required, please combine clinical symptoms and other testing methods.

**Note**

1. All operations must be strictly carried out in accordance with the instructions;
2. The various components in the reagent kit should be naturally melted, completely mixed, and briefly centrifuged before use;
3. The reaction solution should be stored away from light;
4. Try to avoid the presence of bubbles during the reaction, and cover the tube tightly;
5. Use disposable suction tips, disposable gloves, and specialized work clothes for each area;
6. Sample processing, reagent preparation, and sample addition should be carried out in different areas to avoid cross contamination;
7. After the experiment is completed, treat the workbench and pipette with 10% hypochlorous acid, 75% alcohol, or a UV lamp;
8. All items in the reagent kit should be treated as contaminants and handled in accordance with the "Biosafety Guidelines for Microbial Biomedical Laboratories"