



HotStart Taq DNA Polymerase (Glycerol-Free)

Product Number: PC17GF

Shipping and Storage

-20°C.

Components

Component	PC17GF
HotStart Taq DNA Polymerase (Glycerol-Free) (5U/μl)	50μl
5×HotStart Taq Buffer with Mg ²⁺ *	1ml×2

Note:*Not suitable for freeze-drying.

Description

This product is a product of Taq enzyme processed by special processes. After chemical modification, the polymerase activity of Taq enzyme is inhibited before heating to high temperature, thereby inhibiting non-specific amplification caused by primer non-specific annealing or primer dimerization under low temperature conditions. This product is suitable for high-specificity PCR reactions, Multiplex PCR, high GC content (>60%), secondary structure, and other genome amplification and large-scale genome amplification detection with strong background.

This product does not contain components such as glycerol that affect the freeze-drying process and can be used for the preparation of freeze-drying reaction systems and product design.

Application

1. High specificity PCR reaction;
2. Complex template amplification;
3. Multiplex PCR;
4. Genomic amplification testing;
5. Fluorescence quantitative PCR.

Unit definition

After thermal activation, the amount of enzyme required to catalyze the incorporation of 10nmol dNTPs into acid insoluble substances within 30 minutes at 74°C is one unit.

Quality control

After multiple column purification, only a clear and single target band was visible in SDS-PAGE gel detection. qPCR method detected no residual Escherichia coli DNA and no contamination of nucleic acid endonucleases and exonucleases.

Suggestion

The PCR product amplified using this reagent has a prominent "A" base at the 3' end, which can be directly cloned into a T vector.

Protocol

1. Common reaction systems (50μl)

5×HotStart Taq Buffer with Mg ²⁺ *	10μl
Upstream primer	0.2-1.0μM(Final Conc.)
Downstream primer	0.2-1.0μM(Final Conc.)



Tinzyme Co., Limited

Email: sales@tinzyme.com

Website: www.tinzyme.com

Tel: +86-755-86134126

WhatsApp/Facebook/Twitter: +86-189-22896756

dNTPs (10mM each)	1µl
Template	1-50ng(Plasmid) 10ng-1µg(Genome)
HotStart Taq DNA Polymerase	0.25µl(1.25U)
ddH ₂ O	Up to 50µl

2. Common PCR reaction programs:

2.1. Two step method:

Step	Cycles	Temperature	Time
Pre denaturation	1	94°C	5min
Denaturation	} 30	94°C	5s
Annealing/Extension		68°C	1Kb/60s
Final extension	1	72°C	5min
Heat preservation	1	4°C	Heat preservation

2.2. Three step method:

Step	Cycles	Temperature	Time
Pre denaturation	1	94°C	5min
Denaturation	} 30	94°C	20s
Annealing		50-60°C	20s
Extension		72°C	1Kb/60s
Final extension	1	72°C	5min
Heat preservation	1	4°C	Heat preservation