



B.HotStart DNA Polymerase

Product Number: PC97

Shipping and Storage

-20°C.

Components

Component	PC97
B.HotStart DNA Polymerase(5U/ul)	50µl
5×HotStart Taq Buffer with Mg ²⁺	1ml×2
dNTPs (10mM each)	200µl

Description

B.HotStart DNA Polymerase is a mixture of anti Taq monoclonal antibodies and Taq DNA Polymerase, suitable for Hot Start PCR. Before high-temperature heating, anti Taq monoclonal antibodies bind to Taq DNA polymerase, inhibiting the activity of the polymerase, thereby inhibiting non-specific amplification caused by non-specific annealing of primers or primer dimerization under low-temperature conditions. The anti Taq monoclonal antibody has denatured in the initial DNA denaturation step of the PCR reaction, so there is no need for a hot start step and can be used under conventional PCR reaction conditions. This product is suitable for high-specificity PCR reactions, Multiplex PCR, high GC content (>60%), secondary structure, and other genome amplification and large-scale genome amplification detection with strong background.

Features

With the advantages of high specificity, sensitivity, and strong stability, the product has a wide range of uses.

Application

1. High specificity PCR reaction;
2. Complex template amplification;
3. Multiplex PCR;
4. Genomic amplification testing;
5. Fluorescence quantitative PCR.

Quality control

After multiple column purification, only a clear and single target band was visible in SDS-PAGE gel detection, and there was no contamination of nucleic acid endonucleases or exonucleases detected.

Storage buffer

Tris-HCl (pH8.0), 20 mM; KCl, 100mM; EDTA, 0.1 mM; DTT, 1 mM; Tween20, 0.5%; Nonidet P-40, 0.5%; Glycerol, 50%。

Unit definition

Using activated salmon sperm DNA as a template/primer, the activity of an acidic insoluble substance is defined as one active unit (U) when consuming 10 nmol of whole nucleotide at 74°C for 30 minutes.

Suggestion

The PCR product amplified using this product has a 3-terminal "A" protrusion and can be directly cloned into a T vector.



Tinzyme Co., Limited

Email: sales@tinzyme.com

Website: www.tinzyme.com

Tel: +86-755-86134126

WhatsApp/Facebook/Twitter: +86-189-22896756

Protocol

1. Common reaction systems (50 μ l) :

5 \times HotStart Taq Buffer with Mg ²⁺	10 μ l
Upstream primer	0.2-1.0 μ M(Final Conc.)
Downstream primer	0.2-1.0 μ M(Final Conc.)
dNTPs (10mM each)	1 μ l
Template	1-50ng(Plasmid) 10ng-1 μ g(Genome)
B.HotStart DNA Polymerase(5U/ μ l)	0.25 μ l (1.25U)
ddH ₂ O	Up to50 μ l

2. Common PCR reaction programs:

2.1. Two step amplification program

Step	Temperature	Time	Cycle
Pre denaturation	94 $^{\circ}$ C	1min	1
Denaturation	95 $^{\circ}$ C	20s	35-40
Annealing/Extension	60 $^{\circ}$ C	Adjust according to product length	

2.2. Three step amplification program:

Step	Temperature	Time	Cycle
Pre denaturation	94 $^{\circ}$ C	1min	1
Denaturation	95 $^{\circ}$ C	15s	} 35-40
Annealing	50-60 $^{\circ}$ C	20s	
Extension	72 $^{\circ}$ C	Adjust according to product length	