



B.HotStart Taq DNA Polymerase (for-Lyo)

Product Number: PC97LY

Shipping and Storage

-20°C.

Components

Component	PC97LY
B.HotStart Taq DNA Polymerase (for-Lyo)	50µl
5×HotStart Taq Buffer with Mg ²⁺ (for-Lyo)	1ml

Description

B.HotStart Taq DNA Polymerase is a mixture of anti Taq monoclonal antibodies and Taq DNA Polymerase, suitable for Hot Start PCR. Before high-temperature heating, anti Taq monoclonal antibodies bind to Taq DNA polymerase, inhibiting the activity of the polymerase, thereby inhibiting non-specific amplification caused by non-specific annealing of primers or primer dimerization under low-temperature conditions. The anti Taq monoclonal antibody has denatured in the initial DNA denaturation step of the PCR reaction, so there is no need for a hot start step and can be used under conventional PCR reaction conditions. This product is suitable for high-specificity PCR reactions, Multiplex PCR, high GC content (>60%), secondary structure, and other genome amplification and large-scale genome amplification detection with strong background.

This product does not contain components such as glycerol that affect the freeze-drying process, and has added freeze-drying protectants and excipients. It can be directly used for the preparation of freeze-drying samples to improve their stability and usability.

Features

1. This product can be used to make freeze-dried powder or freeze-dried balls;
2. High specificity PCR reaction;
3. Complex template amplification;
4. Multiplex PCR;
5. Genomic amplification testing;
6. Fluorescence quantitative PCR.

Unit definition

Using activated salmon sperm DNA as a template/primer, the activity of an acidic insoluble substance is defined as one active unit (U) when consuming 10 nmol of whole nucleotide at 74°C for 30 minutes.

Quality control

After multiple column purification, only a clear and single target band was visible in SDS-PAGE gel detection. qPCR method detected no residual Escherichia coli DNA and no contamination of nucleic acid endonucleases and exonucleases.

Suggestion

The PCR product amplified using this reagent has a prominent "A" base at the 3' end, which can be directly cloned into a T vector.

Protocol

1. Sample preparation and use
 - 1.1. Prepare freeze-dried samples of B.HotStart Taq DNA Polymerase (for-Lyo) and 5×HotStart Taq Buffer with Mg²⁺

(for-Lyo) as needed;

1.2. Dissolve the prepared sample in RNase Free Water before use. The remaining reagents after dissolution can be stored at -20°C, and can be used directly after melting in the next experiment.

2. Common reaction systems (50µl) :

5×HotStart Taq Buffer with Mg ²⁺ (for-Lyo)	10µl
Upstream primer (10µM)	0.2-1.0µM(Final Conc.)
Downstream primer (10µM)	0.2-1.0µM(Final Conc.)
dNTPs (10mM each)	1.0µl
Template	1-50ng(Plasmid) 10ng-1µg(Genome)
B.HotStart Taq DNA Polymerase (for-Lyo)	0.25µl (1.25U)
ddH ₂ O	Up to 50µl

3. Common PCR amplification programs:

3.1. Two step amplification program:

Cycle	Temperature	Time
1	94°C	1min
35-40	95°C	20s
	60°C	Adjust according to product length

3.2. Three step PCR amplification program:

Cycle	Temperature	Time
1	94°C	1min
35-40	95°C	15s
	50-60°C	20s
	72°C	Adjust according to product length

Note

Both B.HotStart Taq DNA Polymerase (for-Lyo) and 5×HotStart Taq Buffer with Mg²⁺ (for-Lyo) in this product already contain freeze-dried excipients. If other freeze-dried excipients need to be added, they need to be evaluated based on actual test data before adding.