



Super Start Probe qPCR Mix (UDG, Lyo)

Product Number: PCK006LY

Shipping and Storage

-20°C.

Components

Component	PCK006LY
2×Super Start Probe qPCR Mix (UDG, Lyo)	5 vials
ROX I	200µl
ROX II	200µl
RNase Free Water	1ml×10

1. The Real Time PCR instrument calibrated with ROX I includes:
ABI 5700/7000/7300/7700/7900/7900 HT/7900 HT Fast, ABI StepOne/StepOnePlus and other instruments.
2. The Real Time PCR instrument calibrated with ROX II includes:
ABI 7500/7500 Fast; ABI ViiA7; ABI Q6, ABI Quant Studio 6/7 Flex; Stratagene MX4000/MX3005P/MX3000P and other instruments.
3. Real Time PCR instruments that do not require ROX calibration include:
Bio-Rad CFX96/CFX 384/iCycler iQ5/iQ/My iQ5/MiniOpticon/Opticon/Opticon 2/Chromo4 ; Eppendorf Mastercycler ep realplex /realplex 2s; Cepheid SmartCycler; Illumina Eco qPCR; Roche Applied Science LightCycler 480; Thermo Scientific PikoReal Cycler; Qiagen/Corbett Rotor-Gene Q/Rotor-Gene 3000/ Rotor-Gene 6000 and other instruments.

Description

Super Start Probe qPCR Mix (UDG, Lyo) is a reagent that uses TaqMan probe method for Real Time PCR; This product is a freeze-dried powder made by pre mixing HotStart Taq DNA polymerase, UDG enzyme, optimized reaction buffer, dNTPs, protectants, and excipients. The freeze-dried powder has a uniform appearance, a block like or sponge like structure, and dense pores. Simply add RNase Free Water to quickly dissolve into a 2×concentration premixed reagent. During the experiment, primers, probes, templates, and water are added to prepare the PCR reaction solution.

The introduction of dUTP/UDG anti pollution system into the formulation can effectively remove pollutants present in the system. At the same time, when the reaction system is heated to 90-95 °C, UDG enzymes are completely inactivated, which can maintain cDNA integrity.

This product is prepared using a special process to produce HotStart Taq DNA polymerase, combined with a specially optimized reaction and freeze-drying formula to produce a freeze-drying reagent. It can effectively inhibit non-specific reactions, improve reaction specificity and amplification efficiency, and accurately quantify over a wider range. And equipped with ROX Reference Dyes of different concentrations, used for correcting inter well fluorescence signal errors with different instruments.

Features

1. Adding RNase Free Water can quickly dissolve;
2. No need for cold chain transportation, can be stored at room temperature and 2-8°C;
3. High specificity: The HotStart Taq DNA polymerase prepared by a special process greatly improves the specificity of PCR amplification;
4. Efficient: The carefully formulated Real Time PCR specific 2×SuperMix has higher amplification efficiency and sensitivity;
5. Quick: The necessary reagents for PCR reaction are collected in one tube, and the reaction system can be prepared in a few minutes.



Quality control

Purity testing: All components have been tested and found to have no residual endonucleases, exonucleases, or nucleic acid residues. The moisture content of freeze-dried powder is less than 3%.

Suggestion

1. Please use new (non polluting) gun heads, microtubes, etc. for the preparation and packaging of the reaction solution to avoid contamination.
2. Primer design: Generally, in order to increase sensitivity and amplification efficiency while suppressing non-specific reactions, the shorter the amplification target sequence, the better, usually below 200bp.

Protocol

1. Dissolve freeze-dried products

Add 1ml of RNase Free Water to each bottle to dissolve the freeze-dried product. The dissolved product is 2×Super Start Probe qPCR Mix (UDG, Lyo); The dissolved reagent can be stored at -20°C and can be used directly after melting in the next experiment.

2. Common reaction systems (20µl)

2×Super Start Probe qPCR Mix (UDG, Lyo)	10µl
Upstream primers	0.2-1.0µM(Final Conc.)
Downstream primers	0.2-1.0µM(Final Conc.)
Probe (10µM)	0.4µl
Template	1-2µl
ROX (to be selected based on machine model)	0.4µl
RNase Free Water	Up to 20µl

3. Common PCR cycles*

*The two-step method can achieve high specificity, while the three-step method can achieve high amplification rate.

3.1. Two step amplification program:

Temperature	Time	Cycles
50°C	2min	
95°C	5min	
95°C	10s	} 45-50
60°C	40s	

3.2. Three step PCR amplification program:

Temperature	Time	Cycles
50°C	2min	
95°C	5min	
95°C	10s	} 45-50
50-60°C	40s	
72°C	30s	