

## Recombinant PNGase F

**Product Number: PSE109**

---

### Shipping and Storage

1. Storage temperature:  $-20\pm 5^{\circ}\text{C}$  for storage
2. Validity period: 2 years.

### Component

Component	PSE109
PNGase F	30 $\mu\text{L}$
10 $\times$ Buffer 1	1mL
10 $\times$ Buffer 2	1mL
NP-40	1mL

### Description

N-glycosidase F (PNGase F) is an effective amidase produced by recombinant Escherichia coli with cloned genes encoding the enzyme. It can cleave the N-sugar chain on glycoproteins, helping to generate carbohydrate free peptides and oligosaccharides with di-N-acetyl chitosan units. N-glycosidase F is the most effective method for removing almost all N-sugar chains from glycoproteins. The product can remove all complex, heterozygous, and high mannose sugar chains from antibodies and their fusion proteins, but if there are alpha 1,3-linked fucose residues in the core structure (often expressed in plant and insect cells as immunoglobulins), they cannot be cleaved. Obtaining accurate N-sugar chain distribution in the shortest possible time is crucial for effective process control. N-glycosidase F is an optimized reagent that can rapidly deglycosylate antibodies and fusion proteins within a few minutes. All N-sugar chains can be released quickly and without preference, and can be directly subjected to downstream chromatography or mass spectrometry analysis to quickly determine the glycosylation of antibodies.

### Active definition

Unit enzyme activity is defined as the amount of enzyme required to remove over 95% of sugar chains from 10 $\mu\text{g}$  denatured RNase B within 1 hour in a 10 $\mu\text{L}$  reaction system at 37 $^{\circ}\text{C}$ .

### Protocol

#### 1. Reaction conditions under denaturing conditions:

For some glycosylated antibodies, it is necessary to denature the sample before complete deglycosylation. The steps are as follows:

- 1.1. Dissolve 1-20 $\mu\text{g}$  of antibodies or glycoproteins in deionized water, dilute to 9 $\mu\text{L}$  with deionized water, and then add 1 $\mu\text{L}$  of 10  $\times$  Buffer 2, resulting in a total volume of 10 $\mu\text{L}$ .
- 1.2. Incubate at 100 $^{\circ}\text{C}$  for 10 minutes, then cool on ice and centrifuge for 10 seconds.
- 1.3. Add 2 $\mu\text{L}$  of 10 $\times$ Buffer 1, 2 $\mu\text{L}$  of 10% NP-40, and 6 $\mu\text{L}$  of deionized water, resulting in a final volume of 20 $\mu\text{L}$ .
- 1.4. Add 1 $\mu\text{L}$  of PNGase F, gently blow and mix well.
- 1.5. Incubate at 37 $^{\circ}\text{C}$  for 1 hour.
- 1.6. The simplest method to evaluate the degree of deglycosylation is to observe the migration changes of the target protein through SDS-PAGE gel electrophoresis.

#### 2. Reaction conditions under non denaturing conditions:

- 2.1. Dissolve 1-20 $\mu\text{g}$  of antibodies or glycoproteins in deionized water, add 2 $\mu\text{L}$  of 10 $\times$ Buffer 1, and the final volume is 20 $\mu\text{L}$ .
- 2.2. Add 2-5 $\mu\text{L}$  of PNGase F, gently blow and mix well;
- 2.3. Incubate at 37 $^{\circ}\text{C}$  for 4-24 hours.



## **MEBEP TECH(HK) Co., Limited**

Email: [sales@mebep.com](mailto:sales@mebep.com) Website: [www.mebep.com](http://www.mebep.com)

Tel: +86-755-86134126 WhatsApp/Facebook/Twitter: +86-189-22896756

Note: Under denaturing conditions, most substrates can undergo better deglycosylation. Under non denaturing conditions, it may be necessary to increase the amount of PNGase F and prolong the incubation time.

2.4. The simplest method to estimate the degree of glycosylation is to observe the migration changes of the target protein through SDS-PAGE gel electrophoresis.