



Malt Extract Agar

Product Number: MEA317-250g

Product details

Formula of culture medium (per liter); pH7.0 ± 0.1(25°C)

Component	Weight
Malt Extract	30.0g
Soya Peptone	3.0g
Agar	15.0g

Instructions

Weigh 48.0g of this product into 1L of distilled water or deionized water, heat and boil until completely dissolved, divide and sterilize at 121°C for 15 minutes. After sterilization, shake well to prevent the agar from settling at the bottom of the vessel and solidifying. Set aside for later use.

Application

Used for the cultivation of mold in the evaluation of the anti-corrosion effect of plastic preservatives.