



Metal bead 2.8mm

Product Number: MTB-201-2.8

Storage

Store at room temperature.

Components

Product Number	Product Name	Package
MTB-201-2.8	Metal bead 2.8mm	5500 beads per pack
—	Manual	1 copy

Description

MEBEP Metal bead 2.8mm (Stainless Steel Beads for Tissue Homogenizer, 2.8mm) can be used with high-throughput tissue grinders or other devices that grind through oscillation, for grinding and mixing samples such as tissues.

These Metal bead have a diameter of approximately 2.8mm, with approximately 5500 beads per pack.

These Metal bead can be reused 5-10 times. After use, wash them with water, sterilize with high pressure or soak in 70% alcohol, and dry for future use.

To avoid cross-contamination between different samples, it is more recommended to use them once.

Notes:

Although these Metal bead are made of stainless steel, they may rust after a certain number of uses due to high-speed grinding in lysis solution. It is not recommended to continue reusing them at this point.

This product is intended for scientific research use by professionals only and should not be used for clinical diagnosis or treatment, food or pharmaceuticals, and should not be stored in ordinary residences.

For your safety and health, please wear a lab coat and disposable gloves when operating.