



Tinzyme Co., Limited

Email: sales@tinzyme.com

Website: www.tinzyme.com

Tel: +86-755-86134126

WhatsApp/Facebook/Twitter: +86-189-22896756

Protease Inhibitor Cocktail (100×)

Product Number: PIC83

Shipping and Storage

-20°C

Components

Component	PIC83
	1ml
Protease Inhibitor Cocktail (100×)	1ml

Notes

1. Protease inhibitors severely damage respiratory mucosa, eyes, and skin, posing a fatal risk if inhaled, swallowed, or absorbed through the skin.
2. Please wear experimental clothing and disposable gloves during operation.

Protocol

Take out the protein extraction reagent required for the experiment and pre cool it. Add Protein Inhibitor Cocktail in a 1:99 ratio(For example, adding 10μl Protein Inhibitor Cocktail to 990μl extraction reagent),make it 1 in the extraction reagent × Working fluid