

# Tinzyme Co., Limited

*Email: sales@tinzyme.com* Website: www.tinzyme.com Tel: +86-755-86134126 WhatsApp/Facebook/Twitter: +86-189-22896756

# **Magnetic Microspheres**

## **Product Number: MB05HC**

## **Product Description**

MB05HC are monodisperse ferrimagnetic beads, ~400nm in diameter. The beads are composed of silica and magnetic material. The magnetic materials are coated with silica and presenting a bead surface with optimized silica-like chemistry. The magnetic beads ensure rapid magnetic mobility and efficient isolation of nucleic acids (DNA/RNA). The beads also feature a low sedimentation rate and favorable reaction kinetics, making them particularly suited for automated assays.

#### **Description of Materials**

- 1. Diameter: ~400nm
- 2. Structure: SiO<sub>2</sub>/Fe<sub>3</sub>O<sub>4</sub>
- 3. Surface group: -OH
- 4. Additives: Trace surfactants
- 5. Bead Color: Brownish Black


# **Storage and Stability**

Store at 2°C to 25°C for three years.

#### **Recommended application**

Viral nucleic acid co-extraction, swabs (pharyngeal swabs, saliva HPV,feces), PCR product purification, cfDNA

#### **Physical and Chemical Properties**



For Research Use Only



Tinzyme Co., Limited

Email: sales@tinzyme.comWebsite: www.tinzyme.comTel: +86-755-86134126WhatsApp/Facebook/Twitter: +86-189-22896756

# **Limited Product Warranty**

Freezing, drying, centrifuging, etc. of the product may cause irreversible agglomeration of the microspheres, resulting in decreased bio-application performance.