



Mueller Hinton Agar(MHA)

Product Number: MHA01

Product details

Formula of culture medium (per liter); pH7.4 ± 0.2(25°C)

Component	Weight
Casamino Acid	17.5g
Beef Extract	2.0g
Soluble Starch	1.5g
Agar	17.0g

Instructions

Weigh 38.0g of this product, add 1000mL of distilled water, stir and heat to boil until completely dissolved, sterilize at 121°C for 15 minutes under high pressure, cool to around 50°C, pour an equal amount and culture based on sterilization culture (4mm high).

Application

Used for agar dilution and diffusion methods in drug susceptibility testing.