



## Potato dextrose agar(PDA)

Product Number: PDA-250g

---

### Product details

Formula of culture medium (per liter); pH5.6  $\pm$  0.2(25°C)

Component	Weight
Infusion from potatoes	11.0g
Dextrose	20.0g
Agar	15.0g

### Instructions

Weigh 46.0g of this product into 1L of distilled water or deionized water, dissolve at a low temperature, divide and sterilize at 115°C for 30 minutes. After sterilization, shake well to prevent the agar from settling at the bottom of the vessel and solidifying. Cool to around 50°C, pour into a petri dish, and set aside.

### Application

Used for counting mold and yeast.